

## Miltex® Surgical Instrument Stain Remover

Powder formula removes stains, tarnish and restores luster to surgical instruments, autoclaves, ultrasonic cleaners, stainless steel countertops, etc. Ideal for use on stainless steel, chrome, nickel, brass, copper and silver.

- Effective on instruments, autoclaves, ultrasonic cleaners or any stainless steel, tile or plated surface.
- Will not harm instrument surfaces
- Easy, fast, convenient to use
- Economical - a little goes a long way



REF	Description
3-740	MILTEX Surgical Instrument Stain Remover, 3 oz. (85 grams), (case of 12).

## Miltex® Nylon & Stainless Steel Brushes

- Nylon bristles for effective cleaning of instrument box locks
- Stainless steel bristles for cleaning serrations



REF	Description
3-1000	MILTEX Instrument Cleaning Brushes 7-1/4" (18.4cm), Nylon Bristles (Set of 3)
3-1001	MILTEX Instrument Cleaning Brushes 7-1/4" (18.4cm), Stainless Steel Bristles (Set of 3)
21-379	Stainless Wire Brush, 6-1/4" (15.9 cm)

## Miltex® Instrument Care Kit System

A great way to try various Miltex products!



REF	Description
3-800	Spray Lube (3-700), Surgical Instrument Cleaner (3-720), Surgical Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon Brush and Stainless Steel Brush.

## Miltex® Instrument Care Training DVD

25-minute Instructional DVD assists practitioners and staff to properly care for instruments and to keep them in top working condition. DVD offers step-by-step guidance on handling, cleaning, rinsing, inspecting, lubricating and sterilizing surgical instruments



REF	Description
3-810	MILTEX Instrument Care DVD

www.miltex.com  
toll-free phone 866-854-8300  
toll-free fax 866-854-8400  
phone 717-840-9335  
fax 717-840-9347

Please call Customer Service today for additional information, an in-office demonstration, to place an order, or to receive other Miltex catalogs.



INSTRUMENT CARE & CLEANING PRODUCTS

Miltex®, Inc.

All Rights Reserved

0708/10/G/8

INST-CARE-BROC



\*+H834INSTCAREBROC07\*

## Miltex® Instrument Care & Cleaning Products

The Miltex® care and cleaning product portfolio is designed to provide optimum cleaning performance and will help to prolong the service life of valuable surgical and dental instrumentation.

## Miltex® Instrument Prep Enzyme Foam

Ready to use foaming spray for pre-cleaning of soiled instruments and scopes (rigid or flexible). Ideal for applications where instrument cleaning is not immediately available. Simply spray on instruments until ready for decontamination and cleaning.



REF	Description
3-760	Instrument Prep Enzyme Foam, 24 oz. (0.71 liter) spray bottles (case of 12)

(Contact Distributor for availability of eaches)

## Miltex® Spray Lube (Surgical Instrument Lubricant)

Pre-mixed, non-silicone and anti-corrosive formula in a ready to use 8 oz. spray pump bottle. Steam permeable and water soluble. Helps prevent spotting, staining and rusting. Lubricates instrument box locks, screw locks, scissor blades and any other metal on metal moving part.



REF	Description
3-700	Spray Lube, 8 fl. oz. (0.24 liter) bottles (case of 12).

(Contact Distributor for availability of eaches)

## Miltex® EZ-Zyme® All Purpose Enzyme Cleaner

Ideal for cleaning hard to remove proteinaceous (mucous & tissue) and lipid (fatty) materials from the surfaces of surgical and dental instruments, linen, water submersible equipment, or where mechanical scrubbing is impractical.

- Easy to Use
- Neutral pH
- Multi-use as an instrument presoak, ultrasonic cleaner, laundry presoak and evacuation system cleaner



REF	Description
3-750	EZ-Zyme All Purpose Enzyme Cleaner, 3/4 oz. (22 ml) packets (dispenser box of 32). Contents of 1-box makes 32 gallons of solution.
3-755	EZ-Zyme All Purpose Enzyme Cleaner, 1-gallon bottle (3.79 liter). Contents of 1-gallon container makes 170 gallons of solution.

(Contact Distributor for availability of eaches)

## Miltex® Surgical Instrument Cleaner

Super concentrate formula for manual and ultrasonic cleaning or soaking of surgical instruments, accessories and glassware.

- Neutral pH
- Biodegradable
- Rinses free leaves no residue
- Phosphate-free, never stains instruments



REF	Description
3-720	Surgical Instrument Cleaner, 8 fl. oz. (0.24 liter) bottles (case of 12). Contents of one 8 oz. bottle makes 32 gallons of solution.
3-725	Surgical Instrument Cleaner, 1-gallon bottle (3.79 liter). Contents of 1-gallon container makes 512 gallons of solution.
3-726	Surgical Instrument Cleaner, 1/4 fl. oz. (7.4 ml) packets (dispenser box of 80). Contents of 1-box makes 80 gallons of solution.

(Contact Distributor for availability of eaches)

## Miltex® Surgical Instrument Cleaner & Lubricant

A dual-purpose instrument cleaner and lubricant. Chemically engineered for use as a presoak or for manual cleaning, washer/decontaminators, ultrasonic cleaners and cart washers.

- Cleaner & Brighter Instruments
- Removes Stains, Blood and Mineral Deposits
- Biodegradable
- Concentrate, Lower Cost Per Use



REF	Description
3-710	1 Gallon (128 fl. oz.) 3.79L (see label for concentration rates)

## Tip-It® Instrument Guards

Miltex® Tip-It® Instrument Guards protect instrument tips and edges from damage during handling and sterilization.

- PVC medical grade plastic
- Latex-free material
- Large holes (vented styles) to reduce vacuum and moisture build up for easy removal from instrument



REF	REF	REF	REF	Dimension	Size Code	Color
Clear Vented	Tinted Vented	Solid Vented	Solid Non-Vented			
3-2501C	3-2501T	3-2501V	3-2501	1/16" X 3/4" (1.6mm X 19mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" X 3/4" (2mm X 19mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" X 3/4" (2.8mm X 19mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" X 1" (3.2mm X 25.4mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" X 1" (4.8mm X 25.4mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" X 3/4" (10.3mm X 19mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" X 3/8" X 1" (1.6mm X 9.5mm X 25.4mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" X 5/8" X 1" (2mm X 16mm X 25.4mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" X 1" X 1" (3.2mm X 25.4mm X 25.4mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5		

\* Color does not apply to clear vented pattern